

Applied Biosystems

7300 Real-Time PCR System

A *Real* Affordable Approach to Real-Time PCR

- Four-color detection provides the flexibility to perform a variety of applications, including gene expression analysis, pathogen quantitation, SNP genotyping, isothermal, and plus/minus assays that utilize internal positive controls
- Powerful, versatile software includes plate set-up wizards that guide you through experimental set-up; advanced data-viewing capabilities including data sorting and custom graphing options, simple export to PowerPoint, and automated analysis tools that make data processing simple and straightforward
- Precision optics and a chargecoupled device (CCD) camera, combined with a sophisticated multicomponenting algorithm, provide highly accurate, reproducible, and reliable results
- Latest generation, Peltier-based, thermal cycling system accommodates both standard 96-well plates and 0.2 mL tubes



The Applied Biosystems Advantage

Built on over 10 years of real-time expertise, the Applied Biosystems 7300 Real-Time PCR System is an affordable platform for the detection and quantification of nucleic acid sequences that will not compromise your data quality or dye choice flexibility. The 7300 Real-Time PCR System combines thermal cycling, fluorescence detection, and application-specific software to measure the cycle-by-cycle accumulation of PCR products in a single-tube, homogeneous reaction.

Quantitative results are available immediately upon completion of PCR, with no need to run gels, purify PCR products, or perform any post-PCR manipulations. Real-time PCR runs are completed in under two hours, using 96-well plates and tubes (individual or 8-strip), with a supported volume range of 20 – 100 μ L. Compared with manual PCR quantitation techniques such as Northern blotting or RNase protection assays, real-time PCR offers enormous time-savings, greater sensitivity, superior precision, and a larger dynamic range. This high-quality platform from the leader in real-time PCR systems provides data you can trust at a price you can afford.

Real-Time PCR Applications

The 7300 system supports many realtime quantitative PCR applications including gene expression analysis using relative quantitation (RQ) assays, and absolute quantitation using standard curves. In addition, the system allows for qualitative, post-PCR detection of nucleic acids for allelic discrimination (SNP genotyping) assays and plus/minus assays that use internal positive controls.

Fluorescence Detection

All sample wells in the 7300 system are illuminated with a tungsten-halogen lamp. Fluorescence emission is detected through four filters on to a CCD camera. The emission filters are optimized for use with FAMTM/SYBR[®] Green I, VIC[®]/JOETM, NEDTM, TAMRATM, and ROXTM fluorescent dyes.

Sequence Detection Software

Instrument software for the Applied Biosystems 7300 Real-Time PCR System runs on the Windows XP[®] operating system and provides instrument control, data collection, and data analysis. Powerful and user-friendly, sequence detection software includes the following features:

- Plate set-up wizards for easy experimental design, even with complex multicolor assays
- Real-time monitoring of amplification growth curves enables you to view run progress
- Auto-baseline and auto-threshold for simplified data analysis
- Absolute quantitation of nucleic acid targets with the ability to simultaneously analyze multiple standard curves on a single plate
- Optional relative quantitation (RQ) study software with powerful dataviewing capabilities allowing the simultaneous analysis of up to ten 96-well plates containing gene expression data

Instrument Specifications

Instrument	7300 System
Thermal cycling system	Peltier-based, 96-well block
Sample Ramp Rate	+/-1.1°C/sec
Peak Block Ramp Rate	1.6°C/sec
Temperature Range	4°C – 100°C
Temperature Accuracy	±0.25°C of setpoint/display temperature, measured at 3 minutes after clock start.
Temperature Uniformity	±0.50°C, 30 seconds after clock start
Optical system	Single excitation, four emission filters, and CCD camera
Calibrated Dyes at Installation	SYBR [®] Green I, FAM TM , VIC [®] , JOE TM , NED TM , TAMRA TM , ROX TM
Passive Reference Dyes	ROX TM or any calibrated dye.
Data Collection	Data collected in all 4 filters for all wells regardless of plate setup. Plate setup may be altered after run completes.
Quantitative PCR run time	< 1 hour 50 minutes
Supported Volumes	20 – 100 µL
Supported Consumables	<ul style="list-style-type: none">• Standard Optical 96-well plates• 8-strip 0.2mL tubes• 0.2mL tubes• Optical Adhesive Covers• Optical Flat Caps

- Automated SNP genotype calling capability with intuitive graphical output and quality-value assignment
- Simple dissociation curve data collection and viewing
- Tool tips for easy identification of sample wells when viewing amplification curves or SNP genotyping plots
- Lamp-life monitoring and instrument diagnostics provide confidence in your instrument's performance
- Data sorting and filtering for easy viewing and reporting, along with customizable graphs for targets or samples of interest
- Versatile system can be used as a plate reader or a regular thermal cycler if needed
- Simple correction of plate setup mistakes without losing collected data

Computer Specifications

Applied Biosystems supplies a DellTM Business Line computer (notebook or tower) for use with the 7300 system. For the latest computer specifications, please visit the Applied Biosystems Web site at www.appliedbiosystems.com.

Installation Specifications

Using the TaqMan[®] RNase P Instrument Verification Plate, the 7300 system can distinguish between samples containing 5,000 and 10,000 template copies with a confidence level of 99.7%.

Demonstrated Performance

The 7300 system has been demonstrated to achieve the following performance targets:

- 9 logs of linear dynamic range
- Detection of 10 starting copies of a DNA template in a 50 µL reaction for a single reporter TaqMan® assay with a confidence level of 99.7%

Reagents and Disposables

A complete line of reagents including TaqMan® Universal PCR Master Mixes and SYBR® Green I Master Mixes, and disposables including tubes and 96-well plates, are available for use with the 7300 system.

TaqMan® Genomic Assays

Applied Biosystems provides preformulated, ready-to-use, quality-tested, 5' nuclease TaqMan® probe-based assays for use with the 7300 system (see table above right).

Service and Warranty

Purchase of the instrument includes a one-year limited warranty on parts and labor, plus an installation package that includes set-up and calibration of the instrument from our highly trained Service Support team.

Support

Applied Biosystems technical specialists and scientists provide worldwide applications support and service.

TaqMan® Genomic Assay Products

Applied Biosystems provides preformulated, ready-to-use, quality-tested, 5' nuclease TaqMan® assays for use with the 7500 and 7500 Fast systems.

Assay Details	Applications	
Assay Formats	Gene Expression	SNP Genotyping
TaqMan® Genomic Assays	Yes	Yes
Custom TaqMan® Genomic Assays	Yes	Yes
Genome Availability	Number of Assays	
Human	> 204,000	> 1,990,000
Mouse	> 179,000	N/A
Rat	> 128,000	N/A
D. Melanogaster (Drosophila)	> 38,000	N/A
A. Thaliana (Arabidopsis)	> 95,000	N/A
URL (assays can be ordered online only)	www.allgenes.com	www.allsnps.com

Instrument and Computer Dimensions

Dimension	7300 System	Notebook	Tower
Width	34 cm (13.39 in.)	32 cm (12.4 in.)	18 cm (7.1 in.)
Depth	45 cm (17.72 in.)	26 cm (10.1 in.)	45 cm (17.6 in.)
Height	49 cm (19.29 in.)	3 cm (1.2 in. closed)	42 cm (16.7 in.)
Weight	29 kg (64 lb)	2.1 kg (4.7 lb)	32 kg (70 lb)

ORDERING INFORMATION

Description	P/N
7300 Real-Time PCR System with Notebook Computer	4351101
7300 Real-Time PCR System with Tower Computer	4351103
RQ Study Software Upgrade	4350814

NOTICE TO PURCHASER:

The Applied Biosystems 7300 Real-Time PCR System is a real-time thermal cycler covered by one or more of US Patents Nos. 6,814,934, 5,038,852, 5,333,675, 5,656,493, 5,475,610, 5,602,756, 6,703,236, 6,818,437, 7,008,789, 6,982,166, and 6,677,151 and corresponding claims in their non-US counterparts, owned by Applied Biosystems Corporation. No right is conveyed expressly, by implication or by estoppel under any other patent claim, such as claims to apparatus, reagents, kits, or methods such as 5' nuclease methods. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

© 2006 Applied Biosystems. All rights reserved. Applied Biosystems and VIC are registered trademarks and AB (Design), FAM, iScience, iScience (Design), JOE, ROX, and TAMRA are trademarks of Applied Biosystems Corporation or its subsidiaries in the US and/or certain other countries. TaqMan is a registered trademark of Roche Molecular Systems, Inc. Dell is a trademark of Dell Corporation. Microsoft and Windows XP are registered trademarks of Microsoft Corporation. Pentium is a registered trademark of Intel Corporation. SYBR is a registered trademark of Molecular Probes, Inc. All other trademarks are the property of their respective owners.

Printed in the USA, 08/06, Publication 117SP04-05



Headquarters

850 Lincoln Centre Drive | Foster City, CA 94404 USA
Phone 650.638.5800 | Toll Free 800.345.5224
www.appliedbiosystems.com

International Sales

For our office locations please call the division headquarters or refer to our Web site at
www.appliedbiosystems.com/about/offices.cfm