RIDA[®] AllergyScreen[®]

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RIDA® AllergyScreen® – for in-vitro diagnostic

Test principle: the immuno-blot analysis

The prevalence of allergies is increasing worldwide. In Europe alone up to 15 % of the population is thought to be afflicted.

Advances in allergy-specific laboratory tests are now indispensable tools for the diagnosis of allergies. It is necessary to test a patient's blood serum whenever an allergic response or severe reaction is suspected. Patients show often a wide range of varying symptoms and sensitization patterns against several allergens.

The RIDA® AllergyScreen® system is an in vitro system for proving the presence of specific IgE antibodies in serum. We would like to give you some information about this new screening concept: RIDA® AllergyScreen®.

Allergens, prepared only for in-vitro diagnostic purpose, are bound as test lines to the surface of nitrocellulose membranes which are glued after production in reaction troughs. The RIDA® AllergyScreen® system supports up to 20 different allergens on one membrane.

For the immunological reaction 250 µl of patient serum are pipetted into the trough and incubated at room temperature. During this time, the allergen-specific IgE-antibodies react with the allergen bands and bind to the nitrocellulose membranes via the allergen. Non-bound material is removed by washing.

After this, an anti-human IgE antibody coupled with biotin is added and incubated at room temperature. This binds to the bound specific IgE antibodies in the test fields from the first incubation. Non-bound detection antibodies are removed by washing.

Next, a streptavidin is added which is conjugated with alkaline phosphatase and incubated at room temperature. The streptavidin binds to the biotin from the second incubation in the test fields and to the positive control.

Non-bound streptavidin conjugate is removed by washing. After adding the substrate (BCIP/NBT) and incubating at room temperature,

RIDA® X-Screen (single)

RIDA® quadro-Screen

RIDA® maXi-Screen 2

a specific enzymatic colour reaction of the alkaline phosphatase takes place. Positive bands become visible as a result of the formation of precipitates on the test membrane. The colouration is directly proportional to the specific antibody content of the serum sample.

Using 2 membrane panel holders, able to carry 10 panels each, 20 patient sera can be manually tested in one single test run.

For the performing of the test a shaker with reciprocal or rocking motion is necessary for an even distribution of the reagents over the test membrane.

Evaluation is carried out after complete drying of the test membrane with an instrument including a CCD camera (RIDA® X-Screen, RIDA® quadro-Screen or RIDA® maXi-Screen 2), which takes a photo of the membrane. The colour intensity of each allergen band is evaluated in accordance with the standard curve which is integrated in the software. The specific IgE concentrations are calculated in IU/ml and assigned to EAST classes (0-6).

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The product

Available allergens

There are different panels with specific compositions of up to 20 allergens available. The advantage of RIDA® AllergyScreen® compared to a single allergen test system is the simple analysis of a whole range of allergens in one operation, whereby only 250 µl of serum are necessary for 20 allergens. The assay time is only 2.5 hours.

Panel 1 is a mixture of inhalative and food allergens for an efficient screening of a possible sensitization.

Panel 2 includes the most important inhalative allergens which covers the most important inhalative sensitizations.

Panel 3 contains only food allergens, which are characterized as most common allergens.

Panel 4 contains a selection of 20 inhalative and food allergens taking special consideration of children and their food.

The results are printed out in IU/ml (range < 0,35 - 100 IU/ml) as well as in EAST classes (0-6). In this way, an efficient and economic determination of specific IgE is given compared to a single allergen test system.

Please ask for specific allergen panels for your country !

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Bestell-Nr. A2342

Bestell-Nr. A2442

Composition of AllergyScreen standard panels

Panel 1 inhalative- und food allergens	Panel 2 inhalative allergens	Panel 3 food-allergens	Panel 4 paediatric panel
Pos. control	Pos. control	Pos. control	Pos. control
1 Der. pteronyssinus	1 Der. pteronyssinus	1 Hazelnut	1 Der. pteronyssinus
2 Der. farinae	2 Der. farinae	2 Peanut	2 Der. farinae
3 Alder pollens	3 Alder pollens	3 Walnut	3 Birch pollens
4 Birch pollens	4 Birch pollens	4 Almond	4 Grass pollens
5 Hazel pollens	5 Hazel pollens	5 Milk	5 Cat
6 Grass pollens	6 Oak pollens	6 Egg white	6 Dog
7 Rye pollens	7 Grass pollens	7 Egg yolk	7 Alternaria alternata
8 Mugwort	8 Rye pollens	8 Casein	8 Milk
9 Plantain	9 Mugwort	9 Potatoe	9 α-Lactoalbumin
10 Cat	10 Plantain	10 Cellery	10 β-Lactoglobulin
11 Horse	11 Cat	11 Carrot	11 Casein
12 Dog	12 Horse	12 Tomato	12 Egg white
13 Alternaria alternata	13 Dog	13 Cod fish	13 Egg yolk
14 Egg white	14 Guinea pig	14 Crab	14 Bovine serum albumin
15 Milk	15 Golden hamster	15 Orange	15 Soya bean
16 Peanut	16 Rabbit	16 Apple	16 Carrot
17 Hazelnut	17 Pen. notatum	17 Wheat flour	17 Potatoe
18 Carrot	18 Cladosp. herbarum	18 Rye flour	18 Wheat flour
19 Wheat flour	19 Asp. fumigatus	19 Sesame seed	19 Hazelnut
20 Soya bean	20 Alternaria alternata	20 Soya bean	20 Peanut

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
Analysis

The basis of development of the software is the digital photographic evaluation of Western Blot bands. The RIDA® X-Screen, RIDA® quadro-Screen and RIDA® maXi-Screen 2 take photographs of the test membranes with a CCD camera and the software evaluates the colour intensity of the allergen bands using a mathematical algorithm.

This is achieved by comparing the calculated surface integral of each band with an internal standard curve and grouping these results into

classes (0-6). Classes are directly related to the specific IgE content in IU/ml.

After measurement a printout provides the user with a photo of the membrane, the densitometer curve of the membrane, the concentration data for each allergen band in IU/ml and the according classes. The data are stored and documented in the system, specifically for each patient, and are retrievable. Additional the analysis can be done visually by the use of a colour card.

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Automation

With the help of an immunoblot incubator like the Bee Blot from Bee Robotics 36 RIDA® AllergyScreen® panels can be developed automatically. The combination of Bee Blot with RIDA® maXi-Screen 2 avoids handling failures and the mix-up of panels due to manual transfers. RIDA® maXi-Screen 2 and RIDA® quadro-Screen support online connection to LIS systems offering work list management and transfer of results back to the LIS.

Therefore the RIDA® AllergyScreen® system will become useful for high throughput laboratories leading to enhanced efficiency and ease of performance of allergy testing.

Only 250 µl patients sera have to be pipetted manually into the reaction troughs. After the incubation steps and drying the whole comb with 36 RIDA® AllergyScreen® panels can be transferred into the RIDA® maXi-Screen 2 for automated evaluation.

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R-Biopharm AG was founded in 1988 as a subsidiary of the company Röhm AG in Darmstadt. In 1991 R-Biopharm was acquired by Dr. Ralf Dreher through a management buy out. The main activity is the development, production and marketing of two different product lines.

1. Enzyme immunoassays for the detection of residues in food and feed

This product line contains test systems for the rapid screening for mycotoxins, hormones and anabolics, antibiotics, vitamins, food allergens and some other parameters.

2. Enzyme immunoassays for the detection of human infectious diseases

The product line contains test systems in the field of stool diagnostics, infectious disease, serology and allergy testing.

All products, agri/food diagnostics as well as the clinical diagnostics, are distributed world-wide under the trademarks RIDASCREEN® and RIDA® via an international distributor network or through own daughter companies like in the US (R-Biopharm Inc.), in the UK (R-Biopharm Rhône Ltd), in France (R-Biopharm France), in Italy (R-Biopharm Italia), in Argentina


(R-Biopharm Latinoamérica), in Brazil (R-Biopharm Brasil), in Spain (R-Biopharm España), in China (R-Biopharm Analysis Systems Trading Co. Ltd.) and in Australia (R-Biopharm Australia).

Since November 1996 the quality management of R-Biopharm is certified according to DIN EN ISO 9001 and DIN EN 46001 (medical products). Certification according to the Quality Management Standard ISO 13485 (Medical product) took place in 2003. Since January 1st 2000 R-Biopharm is the world-wide exclusive distributor for the Boehringer Mannheim Roche enzymatic tests.

Your local distributor:

R-Biopharm AG
An der neuen Bergstraße 17, 64297 Darmstadt, Germany
Phone: + 49 (0) 61 51 - 81 02-0
Fax: + 49 (0) 61 51 - 81 02-20 / -40
www.r-biopharm.com, E-mail: info@r-biopharm.de

r-biopharm

RIDA[®] AllergyScreen

Article no: A2142 Panel 1 / ASAN / HVEN / IND /
IR / KO / KZ / MA / NP / OUSA / PE / TW / UY / VE / VIET /

Article no: A2242 Panel 2 / ASAN / CA / DOHA / H / IND /
KO / KSA / KZ / MA / ME / MENA / MOFID / TR / VE /

Article no: A2342 Panel 3 / CA / DOHA / H / IND / KO / KSA / KZ /
MA / ME / MENA / MOFID / RO / SK / TR / VE /

Article no: A2442 Panel 4 / KO / KZ / MA / VE

Article no: A3054

R-Biopharm AG, An der neuen Bergstraße 17, D-64297 Darmstadt,
Germany, Tel: +49 (0) 61 51 81 02-0 / Fax: +49 (0) 61 51 81 02-20

1. Intended use

For *in vitro* diagnostic use. This is an enzyme immunoassay on a nitrocellulose membrane (immunoblot) for the semi-quantitative determination of specific IgE antibodies against a panel of individual allergens in human serum.

2. Summary and explanation of the test

The purpose of the immune system is to defend the body against pathogenic bacteria, viruses and other microorganisms. The defence reaction serves to protect the organism on initial contact with the pathogens and immunise it on repeated contact. All allergic reactions are preceded by an initial contact phase without symptoms during which Class E specific antibodies (IgE antibodies) are formed. On repeated contact with the allergens which trigger the reaction, these IgE antibodies react with the allergens and lead to the release of mediators (from mast cells or mastocytes) such as histamine, leucotrien and prostaglandine etc. which lead to the symptoms of the allergy. When there is an allergic reaction, the allergens causing the reaction can be identified by determining the specific IgE antibodies in the serum. This method can also be used to determine sensitisations without symptoms.

3. Test principle

This test is based on the principle of an enzyme immunoassay on a nitrocellulose membrane (immunoblot). The allergens which correspond to the panel composition are applied to the surface of nitrocellulose membranes. Allergen-specific IgE antibodies which are present in patient samples react with the antigens and attach themselves to biotin-coupled anti-human IgE antibodies during a second step. During the third incubation step, the Biotin attaches itself to a streptavidin conjugated with alkaline phosphatase (Conjugate). The enzyme converts the colourless substrate (BCIP/NBT) to a blue/lilac end product in the last incubation step. Between the single incubation steps a washing step must be performed. The intensity of the blue colour is proportional to the quantity of allergen-specific antibodies in the serum. The evaluation is carried out by using an analysis template or with the RIDA[®] X-Screen / RIDA[®] quadro-Screen or RIDA[®] maXi-Screen 2.

4. Reagents provided

Table 1: Kits for manual execution:

Membrane	10 pieces	RIDA® AllergyScreen test membrane (nitrocellulose membrane) coated with different allergens on 20 test fields in a reaction trough
Wash	20 ml	Wash buffer, x 25 concentrate, Tris / NaCl
Antibody	4 ml	Detector antibody; anti-human IgE antibodies (goat) conjugated with biotin, ready for use
Conjugate	4 ml	Streptavidin conjugate; streptavidin conjugated with alkaline phosphatase, ready for use
Substrate	4 ml	Substrate; BCIP/NBT (bromochloroindolyl phosphate / Nitro Blue Tetrazolium), ready for use

Table 2: RIDA® AllergyScreen Reagent Kit (A3054)

Wash	2 x 20 ml	Wash buffer, x 25 concentrate, Tris / NaCl
Antibody	2 x 4 ml	Detector antibody; anti-human IgE antibodies conjugated with biotin, ready for use
Conjugate	2 x 4 ml	Streptavidin conjugate; streptavidin conjugated with alkaline phosphatase, ready for use
Substrate	2 x 4 ml	Substrate; BCIP / NBT (bromochloroindolyl phosphate / Nitro Blue Tetrazolium), ready for use

5. Storage instructions

The test membranes must be stored in the plastic packaging under cool, dry, dark conditions. The test kit must be stored at 2-8 °C and can be used after opening until the expiry date printed on the label. As long as the diluted wash buffer is stored at 2-8 °C, it can be used for a maximum of 4 weeks. Microbial contamination must be prevented. After the expiry date has been reached, the quality guarantee is no longer valid.

It is imperative that the conjugate is prevented from contaminating the substrate solution because this will discolour the substrate. The substrate must also be protected from direct light in order to prevent decomposition or discolouration due to auto-oxidation. Once the colour has changed, the substrate must not be used.

6. Additional necessary reagents – and necessary equipment

6.1. Reagents

- Distilled or deionised water

6.2. Accessories

- Vortex mixer
- Measuring cylinder (500 ml)
- Micropipette 250 µl
- 500 ml Laboratory wash bottle
- Test membrane holder for holding 10 test-membranes (optional)
- Incubation box for incubation in the dark (system consisting of test membrane holder and incubation box can be obtained from R-Biopharm)
- Horizontal shaker (optional)
- Hairdryer, standard (optional)
- RIDA® X-Screen measuring instrument including software plus personal computer with USB port (optional)
- RIDA® quadro-Screen measuring instrument including computer and software (optional)
- RIDA® maXi-Screen 2 measuring instrument including computer and software (optional)

7. Precautions for users

For *in vitro* diagnostic use only.

This test must only be carried out by trained laboratory personnel. The guidelines for working in medical laboratories must be followed and the instructions for carrying out the test must be strictly adhered to.

Samples or reagents must not be pipetted by mouth and contact with injured skin or mucous membranes must be prevented. When handling the samples, wear disposable gloves and when the test is finished, wash your hands. Do not smoke, eat or drink in areas where samples or test reagents are being used.

Antibody and wash buffer concentrate contain sodium azide as a preservative. This substance must not be allowed to come into contact with the skin or mucous mem-

brane. Explosive metal azides may be produced on contact with lead or copper pipes.

Conjugate contains methylisothiazolone and bromonitrodexan in subtoxic concentrations as a preservative.

If the outer packaging is defective, the individual components must be examined to make sure they are undamaged before use. Kit components must not be used if their individual packaging is damaged or their containers are not sealed.

You are solely responsible for the proper disposal of all components in the kit after use.

All reagents and materials coming into contact with potentially infectious samples must be treated with suitable disinfectants or autoclaved at 121°C for at least 1 hour.

8. Specimen collection and storage

The test has been developed for testing human serum. After blood collection, the blood should be separated from blood clots as soon as possible in order to prevent haemolysis. The samples must be stored cold or frozen until they are tested. Repeated freezing and thawing of the serum and microbial contamination must be avoided at all costs. Using heat-inactivated, lipaemic, haemolytic, icteric or turbid sera can lead to false results.

Table 3: Sample storage

Undiluted serum	
2-8 °C	-20 °C
1 week	>1 week

9. Test procedure

9.1. General information

All reagents, patients' sera and test membranes must be brought to room temperature before use. The reagents must be thoroughly mixed immediately before use. After use, the kit must be immediately returned to storage at 2-8 °C.

The test membranes cannot be used more than once. The reagents and test membranes must not be used if the packaging is damaged or the containers are not sealed.

Components from kits with different lot numbers must not be combined or exchanged.

Reproducible results strongly depend on keeping to the incubation times and temperatures as well as washing the test membranes uniformly.

The test must not be carried out in direct sunlight. The test membranes must only be held by the handle. Contact with the reaction surface must be avoided – do not touch it. The handle of the reaction trough can be marked with patient data (e.g. laboratory numbers) – with water resistant felt-tip pen.

9.2. Preparing the wash buffer

Make up each bottle of wash buffer concentrate [Wash] to 500 ml with distilled water. Any crystals present in the concentrate must be dissolved beforehand by warming in a water bath at 37°C. Fill a laboratory wash bottle with diluted wash buffer.

9.3. First incubation

Take the test membranes [Membrane] out of the packaging according to the number of tests which are to be carried out. For simplicity, the test membrane holder which can hold 10 test membranes [Membrane] can be used. Wet the test membranes completely with diluted wash buffer (squeezing bottle) and wait until no tiny bubbles are rising up. This is most easily achieved by keeping the test membrane holder horizontal and, after filling the test membranes, carefully rocking the holder back and forth several times. After this, empty the test membranes [Membrane], rinse them again briefly with diluted wash buffer and hit it upside down onto an absorbent tissue. Subsequently fill the test membranes [Membrane] with 250 µl patient serum and incubate them at room temperature (20-25 °C) on the horizontal shaker (100 - 120 rpm) for 45 minutes. Please make sure that the membranes are covered completely by the fluid.

9.4. Washing

Rinse off the test membranes [Membrane] over the sink with the diluted wash buffer (squeezing bottle) for at least 5 seconds. Hold the test membranes [Membrane] vertically downwards to avoid splashing the sera onto neighbouring test membranes. Direct the jet of wash solution over the test membranes several times. After this, fill the [Membrane] with diluted wash buffer, shake it back and forth several times and empty. Finally, hold the test membranes [Membrane] downwards at an angle again and rinse them with the wash bottle for 5 seconds. After this, empty the test membranes and hit them upside down onto an absorbent tissue.

9.5. Second incubation

Add 5 drops (approx. 250 µl) antibody [Antibody] to each test membrane. Please make sure that the membranes are covered completely by the fluid. Incubate the test

membranes **Membrane** on the horizontal shaker (100 - 120 rpm) at room temperature (20-25 °C) for 45 minutes.

9.6. Washing

Washing – see Section 9.4.

9.7. Third incubation

Add 5 drops (approx. 250 µl) conjugate **Conjugate** to each test membrane. Please make sure that the membranes are covered completely by the fluid. Incubate the test membranes **Membrane** on the horizontal shaker (100 - 120 rpm) at room temperature for 20 minutes.

9.8. Washing

Washing – see Section 9.4.

9.9. Fourth incubation

Add 5 drops (approx. 250 µl) substrate **Substrate** to each test membrane. Please make sure that the membranes are covered completely by the fluid. Incubate the test membranes **Membrane** on the horizontal shaker (100 - 120 rpm) at room temperature in the dark for 20 minutes.

After incubation, terminate the colour reaction by briefly rinsing the test membrane membrane with plenty of distilled water or under running water (tap water). Leave the membranes to dry in the air or use a standard hairdryer to accelerate the drying process. The blue/lilac-coloured background on the test membranes disappears on drying. The test membranes in the reaction trough **Membrane** must not be evaluated until they are completely dry.

10. Quality control – indications of reagent expiry

The test has been correctly carried out if the background has fully disappeared and the positive control shows either a strong band or achieves at least EAST Class 2,5 when evaluated in the RIDA® X-Screen, -quadro-Screen or -maxi-Screen 2.

If the values differ from those required or if the reagent is turbid or the substrate has turned blue before adding it to the test membranes, this may be an indication that the reagents have expired.

If the stipulated values are not met, the following points must be checked before repeating the test:

- Expiry date of the reagents used
- Functionality of the equipment being used (e.g. calibration)

- Correct test procedure
- Visual inspection of the kit components for contamination or leaks – a substrate solution which has turned blue must not be used.

If the conditions are still not fulfilled after repeating the test, please contact your local R-Biopharm distributor.

11. Evaluation and interpretation

11.1. Membrane configurations of RIDA® AllergyScreen Panels 1, 2, 3, and 4

	Panel 1	Panel 2	Panel 3	Panel 4
	Positive control	Positive control	Positive control	Positive control
	Derm. pteronyssinus	Derm. pteronyssinus	Hazelnut	Derm. pteronyssinus
	Derm. farinae	Derm. farinae	Peanut	Derm. farinae
	Alder	Alder	Walnut	Birch
	Birch	Birch	Almond	Grass mixture
	Hazel	Hazel	Milk	Cat
	Grass mixture	Oak	Egg white	Dog
	Rye (Pollen)	Grass mixture	Egg yolk	Alternaria alternata
	Mugwort	Rye (Pollen)	Casein	Milk
	Plantain	Mugwort	Potato	α -Lactalbumin
	Cat	Plantain	Celery	β -Lactoglobulin
	Horse	Cat	Carrot	Casein
	Dog	Horse	Tomato	Egg white
	Alternaria alternata	Dog	Cod	Egg yolk
	Protein	Guinea pig	Crab	Bovine serum albumin
	Milk	Hamster	Orange	Soybean
	Peanut	Rabbit	Apple	Carrot
	Hazelnut	Penicillium notatum	Wheat flour	Potato
	Carrot	Cladospor. herbarum	Rye flour	Wheat flour
	Wheat flour	Aspergillus fumigatus	Sesame	Hazelnut
	Soybean	Alternaria alternata	Soybean	Peanut

The membrane configurations for all other country specific panels to which this package leaflet applies are available with R-Biopharm AG as supplement for each panel.

11.2. Findings for the sera

11.2.1. Optical-visual assessment

The colour intensity on the test fields is directly proportional to the quantity of specific IgE antibodies in the serum of the patient for the allergen concerned.

If a band appears in relation to the membrane background, then specific antibodies are present in the serum. If no colouring appears in the reaction field on the membrane, then RIDA® AllergyScreen has not detected any IgE antibodies which are specific to this allergen.


Figure 2: Grey scale for the semi-quantitative evaluation (classes 0 - 4)

Table 4: Relationship between the class determined and the allergen-specific IgE content of the patient serum

Class	Allergen-specific IgE content
0	none found or hardly exists
1	low
2	increased
3	significantly increased
4	High

11.2.2. Quantification using RIDA[®] X-Screen / RIDA[®] quadro-Screen / RIDA[®] maXi-Screen 2

For the quantification, the test membranes are inserted into the membrane holder of one of the above mentioned measuring instruments and measured using the associated software. The IU/ml are calculated automatically from the measured values and assigned to the test classes 0 - 6. The evaluation is based on a standard curve which is stored in the evaluation software.

Please pay attention to the reader specific instructions when using one of the measuring instruments.

It is essential to make sure that the test which corresponds to the allergy panel concerned is called before the measurement.

Using Table 5, the allergen-specific IgE titers can be read off from the concentrations determined in IU/ml or classes.

Table 5: Relationship between the determined IU/ml, classes and allergen-specific IgE contents of the patient serum

IU/ml	EAST-Class	Allergen-specific IgE content
0.00 - 0.34	0 (0.0 - 0.9)	none found or hardly exists
0.35 - 0.69	1 (1.0 - 1.9)	low
0.70 - 3.49	2 (2.0 - 2.9)	increased
3.50 - 17.49	3 (3.0 - 3.9)	significantly increased
17.50 - 49.99	4 (4.0 - 4.9)	high
50.00 - 99.99	5 (5.0 - 5.9)	very high
≥ 100.00	6 (≥ 6.0)	extremely high

11.3. Documentation

After drying of the test membranes and evaluation in the RIDA[®] X-Screen / RIDA[®] quadro-Screen or RIDA[®] maXi-Screen 2 or by optical-visual means has been carried out, the membranes can be removed from the reaction troughs using tweezers and documented in an operating record.

The measuring data (photo of the test membranes and evaluation) are stored on the hard disk in the PC in a default directory. A data sheet can be printed out from a standard printer connected to the PC / measuring instrument for each serum tested.

12. Limitations of the method

The IgE concentrations determined by using this test system make it possible to say something about the degree of sensitisation of the patient to the individual allergens or mixtures of allergens tested.

They cannot be used to derive a relationship between the determined IgE concentration and occurrence of serious clinical symptoms. The results obtained must always be interpreted in combination with the complete clinical picture.

Because of the absence of national and international standards and because of the possible differences in prick-test solutions and allergen extracts which are used for the *in-vitro* tests, it is possible for discrepancies to occur between the results from the *in-vivo* tests and results from the *in-vitro* tests. IgE titres measured immediately after the appearance of anaphylactic reactions may also be negative or too low. If there are discrepancies between the results from the *in-vivo* and *in-vitro* diagnostics, the test should be repeated after 3-4 weeks. Discrepancies which arise should be investigated by an allergologist with a follow-up *in-vivo* test such as a provocation test. Provocation tests can trigger anaphylactic shock.

False positive test results may be produced by cross reactivity of the antigens being tested with other antigens.

13. Performance characteristics

Intra-assay variation: average 4,5 %

Inter-assay variation: average 4,8 %

In order to determine the sensitivity and specificity, 142 sera were tested (in 881 determinations altogether) as part of a clinical study in a quantitative reference *in-vitro* system and compared with the findings of the RIDA[®] Allergy-Screen. 737 skin-prick test results were also compared with the results of the RIDA[®] AllergyScreen-Tests.

Comparison with the IgE reference system:

Sensitivity: 84.3 %

Specificity: 95.0 %

Accuracy 90.6 %

Comparison with the skin prick test:

Sensitivity: 95.1 %

Specificity: 80.2 %

Accuracy 88.3 %

Literature

- Kersten, W., von Wahl, P.G., Lange, C.E., Wenning, J.: Empfehlungen zur in-vitro Diagnostik allergischer Erkrankungen, [Recommendations for the in-vitro diagnosis of allergic disorders] Allergologie, Jg. 23, No. 6, 304-307 (2000)
- von Wahl, P.G., Kersten, W.: Klinische Studie mit einem neuen in-vitro Testsystem, [Clinical study with a new in-vitro test system] Allergo Journal 8, 107 – 111 (1999)
- Kjellman, N.I.M.: Prediction and prevention of atopic allergy, Allergy 37, 463 - 473 (1982)
- Debelic, M.: Die klinische Bedeutung der Bestimmung von Serum-Gesamt IgE und spezifischem IgE für die Diagnostik und Verlaufsbeobachtung allergischer Krankheitsbilder [The clinical importance of determining total serum IgE and specific IgE in the diagnosis and monitoring of allergenic clinical pictures]. Therapiewoche 29, 2280 - 2295 (1979)
- Urbanek, R.: Papier-Radio-Immuno-Sorbent-Test (PRIST) - IgE-Spiegel bei nicht allergischen und allergischen Kindern, [Paper disc Radio Immuno Sorbent Test (PRIST) - IgE levels for non-allergic and allergic children] Monatsschrift Kinderheilkunde 125, 583 - 585 (1977)
- Kjellman, N.I.M., et al.: Serum IgE levels in healthy children quantified by a Sandwich technique (PRIST), Clinical Allergy 6, 51 - 59 (1976)
- Turner, K.J.: IgE globulins and immunity. Medical Journal of Australia 2, 846 (1974)
- Johannson, S.G.O.: Raised levels of a new immunoglobulin class (IgND) in asthma. Lancet II, 951 - 953 (1967)
- Ishizaka, K., Ishizaka, T., Hornbook, M.M.: Physico-chemical properties of human reaginic antibody. IV. Presence of a unique immunoglobulin as a carrier of reaginic activity. Journal of Immunol 97, 75 (1966)
- Bennich, H.H., et al.: Immunoglobulin E, a new class of human immunoglobulin. Bulletin World Health Organisation 38, 151 (1964)