
## SON®

ur requirements



DiaSoria LIAISON°



All reagents are ready to use

for all Eventualities

Wide product menu both for routine

Master curve with 2-point calibration

of Lab Resources

Throughput up to 180 tests/hour

Long walk-away time

Short time to first result

Short hands-on time

and speciality assays

Continuous reloading

On board reagent cooling

Long on board stability

Barcode identification

For effective use

 Ready to use reagent integral Up to 15 assays on board

 15 integrals may be loaded at the same time

### The **LIA SON** Family Collection

HYPERTENSION Direct Renin Aldosterone

SEPSIS BRAHMS PCT®

Anti-HBs Anti-HBc HBc IgM HBeAg Anti-HBe Anti-HAV

HAV IgM HCV Ab

HIV Ab/Ag

Treponema Screen

HTLV I/II Anti HBs II Chagas\*\*

EBV IgM VCA IgG EBNA IgG

The Diagnostic Specialist

VIRAL HEPATITIS AND HIV

HBsAg HBsAg Confirmatory test

#### **Dedicated to Content**





















Borrelia burgdorferi lgG Borrelia burgdorferi lgM









**MEASLES & MUMPS** Measles IgG Measles IgM



Chlamydia T. IgG\* Chlamydia T. IgA\*



. difficile GDH C. difficile Toxin A and B

\* Soon Availab

Product availability subject to required regulatory approval























**STOOL DIAGNOSTICS** 



DiaSorin S.p.A. Via Crescentino 13040 Saluggia (VC) - Italy Tel. +39.0161.487526 Fax: +39.0161.487670 www.diasorin.com E-mail: info@diasorin.it



he LIAISON® analyzer combines flexible operating mode with a broad assay menu and leads to high efficiency due to its long walk-away time. The possibility to run up to 15 different assays at a time provides rational everyday testing even for low numbers of samples. Moreover, a throughput of up to 180 determinations per hour allows LIAISON® to cope easily with a very high volume of samples. It is easy to use due to its touchscreen monitor, clear software structure and the Reagent Integral which contains all assayspecific reagents.



- A stacker holds the reaction modules
- The stacker can hold up to 120 reaction modules or 720 determinations
- Easy reloading even during operation is assured




- Random access, STAT or Batch testing modes
- Rack system for loading up to 144 samples
- All common primary tubes accepted in 3 different racks
- Continuous reloading
- Flexible dilution protocol
- Automatic dilution
- Auto reflex testing
- Barcode identification

# LIAISON®

Meeting all your requirements




LIAISON® work station





LIAISON'





## System Features

#### MODES OF OPERATION

- Random Access
- Batch
- STAT function

#### SAMPLE LOADING

- 144 primary tubes of varying size in racks of 12 samples each
- Continuous reloading

#### REAGENT LOADING

- 15 assays available, allowing up to 1,500 determinations
- Continuous reloading
- Reagent cooling

#### REAGENT CONCEPT

- All assay components are provided ready to use in the Reagent Integral
- Super paramagnetic microparticles
- Chemiluminescence

LIAISON"

#### CALIBRATION CONCEPT

- Master curve procedure with 2-point calibration
- Stable calibration for 2 to 6 weeks

#### PERFORMANCE DATA

- Maximun throughput: 180 tests/hour
- Average throughput: 100-140 tests/hour
- Walk-away time: at least six hours or 720 tests
- Time to first result: 17 mins

#### **FURTHER SYSTEM FEATURES**

- Clot detection
- Barcode concept for samples and reagents
- Continuous reagent and liquid inventory
- Variable dilution protocols (assay- or patient-specific)
- Auto reflex testing
- Automatic dilution of high-titre samples
- Constant incubation temperature of 37 °C
- Operation at ambient temperatures between 15 °C and 30 °C

#### HARDWARE AND SOFTWARE

- IBM-compatible PC
- Colour monitor (touchscreen)

Width/depth/height: 136x66x63 cm (needs bench top)

computer (ASTM protocol)

Bi-directional connection to host

**EDP CONNECTION** 

• Weight: 95 kg

**DIMENSIONS** 



### Reagent integral



- All components, ready to use in one Integral
- Continuous reloading
- Barcode identification



